Machine-learning assisted cavity quench identification at the European XFEL. THPOPA26

LINAC 2022, Oral poster

Branlard, Julien (DESY) Liverpool, 1.9.22





HELMHOLTZ RESEARCH FOR GRAND CHALLENGES


MOTIVATION


Can we lower false positive rates of quench detection at European XFEL?

European XFEL

- 17.5 GeV pulsed electron FEL
- 800x 1.3 GHz SRF cavities
- In operation since 2017
- Hamburg Germany








Data from August 2022


DESY. | Machine-learning assisted quench identification at the European XFEL | Julien Branlard, 1.9.22 | Linac'22

MOTIVATION

Can we lower false positive rates of quench detection at European XFEL?

- Software-based quench detection is successful at catching SRF cavity quenches but also triggers for other faults
- Can we make use of machine learning techniques to help identify REAL and "FAKE" quenches ?
- Make use of data snapshots generated at each trip
 - build up on existing datasets
- Make use of well-known cavity model to compare measured probe and model prediction






1. Introduce a new metric

Residual and Generalized Likelihood Ratio


- Quantifies deviation of the cavity probe from model
- GLR provides very distinct signatures for different kinds of trips



DESY. | Machine-learning assisted quench identification at the European XFEL | Julien Branlard, 1.9.22 | Linac'22

2. Training set on QUENCH – NOT QUENCH


453 trips reviewed by RF expert



DESY. | Machine-learning assisted quench identification at the European XFEL | Julien Branlard, 1.9.22 | Linac'22

3. Run the GLR analysis on all new trip data snapshot

Daily cron job since September 2021



Algorithm accuracy $a = \frac{TP + TN}{TP + TN + FP + FN}$ TP : true positive (accurately detected a quench) TN : true negative (accurately recognized a trip was not a quench) FP : false positive (a "fake" quench) FN : false negative (algorithm failed to identify a real quench)

	TP	TN	FP	FN	а	
QDS	55	56	10	3	~90%	
GLR	55	65	1	3	~97%	
				improved accuracy		

DESY. | Machine-learning assisted quench identification at the European XFEL | Julien Branlard, 1.9.22 | Linac'22

See you at poster THPOPA26

Contact

DESY. Deutsches Elektronen-Synchrotron

www.desy.de

Julien Branlard MSK, DESY julien.branlard@desy.de +49 (0)40 8998 1599

